

MEYRA

iGO



FANTASTIC MANOEUVRABILITY

PRECISE CONTROL
EASY TRANSFERS



Compact and highly manoeuvrable

Great for the indoors—even small apartments



Precise control

Wheelchair controller and max. speed 6kph perfect for safety-critical environments like nursing homes



Easy access, easy transfers

Low open platform at the front, rotating seat for side transfers



Stable and sturdy

User weight up to 120kgs and five wheels—move with confidence, move with safety

Compact and manoeuvrable

THE INDOORS ARE NO OBSTACLE

The Meyra iGo micro is the ideal indoor mobility solution. Thanks to its compact dimensions and tight turning circle, it is very manoeuvrable—critically important for confined indoor spaces. The iGo is lightweight and can be taken apart into three sections for easy transportation. With precise control and five wheels, the Meyra iGo is safe and stable.



Can be disassembled for transport



Tight turning circle



5 wheels for stability and safety



Precise control with wheelchair controller

TECHNICAL DATA

iGO

Length overall	76 cms	Maximum obstacle height	15 mm
Width overall	58 cms	Turning circle	50 cms
Seat width	55 cms	Motor	200 W
Seat depth	40 cms	Batteries	2 x 18ah
Seat height (adjustable)	40 - 55 cms	Driving range	Approx. 17 kms ¹
Backrest height	40 cms	Permitted gradients/ descents	3 degrees/6 degrees
Transport height <i>folded steering column, no seat</i>	47 cms	Empty weight	47 kgs
Front wheel	110 mm	Weight of heaviest piece when disassembled	Approx. 26 kgs
Rear wheel	220 mm	Maximum client weight	120 kgs
Ground clearance	5 cms	Maximum speed	6 kph

¹ distance will be affected by client weight, terrain and battery condition

MEYRA
AUSTRALIA

MEYRA Australia
5 Boulder Rd Malaga
Western Australia 6090



Tel (08) 9248 3560
Fax (08) 9209 1999



sales@meyraaustralia.com.au
www.meyraaustralia.com.au

MEYRA We reserve the right to make technical changes and cosmetic changes to our products. These changes may not be reflected in this brochure.